



FARROW

MEDICAL INNOVATIONS



FarrowWrap Trim-To-Fit™

Better Compression Made Simple™

And Introducing New FarrowWrap™ STRONG Fabric


The Problem: Edema, Lymphedema, Venous Ulceration

- » Long-stretch products are difficult to apply and MAY NOT be the best option
- » Compression stockings and long-stretch products (TED hose, elastic wraps) apply a narrow range of compression to the limb
- » The garment must have enough compression or it will stretch and the limb will swell
- » Garments with enough compression are tight and are often difficult to put on and take off
- » Long-stretch products can hinder blood flow if incorrectly applied
- » Patients struggle donning compression stockings by themselves, even with donning-assist devices

The Solution: FarrowWrap™ Short-Stretch Dynamic Compression



FarrowWrap TTF™



FarrowWrap TTF™ Leg


Free online video inservice for simplified product customization at FarrowMedical.com

Introducing FarrowWrap TTF™ (Trim-To-Fit)


- » Short-stretch technology in an easy to apply ready-made garment
- » Customizable product using our proprietary TTF innovation
- » Available in LITE (20-30mmHg) and STRONG (30-40mmHg) fabrics
- » Simple Velcro® bands adjust as needed without the stress and difficulty of stockings
- » Improved compliance – product is fast, simple to don and adjusts for total comfort

Introducing FarrowWrap Trim-To-Fit™

- » **Can be trimmed to adapt** to extremity as it reduces
- » **Will not fray or roll** when cut
- » **One garment** for active treatment and maintenance
- » **Machine washable** (gentle cycle) and lay flat to dry
- » **Best environmentally friendly** solution for wound care



FarrowWrap TTF™ Leg



Trimming a Footpiece

FarrowWrap™ Technology... Who Benefits?*

- » Patients with weak hand strength
- » Patients with back problems
- » Post-operative patients
- » Patients with large fluctuations in their swelling such as Congestive Heart Failure patients
- » Patients with atypical limb anatomy

* *UT Medical School study shows patients with low Activities of Daily Living (ADL) scores unable to don compression stockings could don the FarrowWrap™ themselves. Reprints available upon request.*

Choose the FarrowWrap™ Product to Fit Your Need

FABRIC

INDICATION

TYPES

LIMB

UPPER EXTREMITY

Arm

Hand Gauntlet

LOWER EXTREMITY

Thigh-High

Knee-High

OPTIONS

FarrowHybrid™ I

– A Class I knee high compression stocking providing 20-30mm compression below level of ankle, with no compression above the ankle. May be used in lieu of a FarrowWrap™ footpiece.

FarrowHybrid™ II

– A Class II knee high compression stocking providing 30-40mm compression below level of ankle, with no compression above the ankle. May be used in lieu of a FarrowWrap™ footpiece.

CLASSIC

Moderate to Heavy
Compression

Custom
Off-the-Shelf

STRONG

Moderate to Heavy
Compression

Trim-to-Fit

LITE

Mild to Moderate
Compression

Custom
Off-the-Shelf
Trim-to-Fit

Call Today for Your Sample Kit!

Free Online Video Inservice at FarrowMedical.com

Farrow Medical Innovations
801 N. Bryan Ave.
Bryan, TX 77803

Toll Free: 877-417-5187 | Fax: 979-775-5202

Pat. Pending
© 2010 FMI. All Rights Reserved.

FarrowWrap Trim-To-Fit™ Better Compression Made Simple™